

Potassium, Urine, Random

CPT Code: 84133

Order Code: C2702

ABN Requirement: No

Specimen: Urine

Volume: 5.0 mL

Minimum Volume: 1.0 mL

Container: Urine (Yellow top, no preservative) tube

Collection:

Urine:

1. Collect and label sample according to standard protocols.
2. Transfer urine from a clean urine cup into a Yellow Top tube using the vacutainer transfer device included with the Yellow Top tube.
3. Cap tube tightly.
4. Label sample according to standard protocols.

Transport: Store urine at 2-8°C after collection and ship the same day per packaging instructions provided with the Cleveland HeartLab shipping box.

Stability:

Ambient (15-25°C): 14 days

Refrigerated (2-8°C): 14 days

Frozen (-20°C): stable

Deep Frozen (-70°C): stable

Causes for Rejection: Specimens other than random urine without preservatives; improper labeling; samples not stored properly; samples older than stability limits

Methodology: Ion Selective Electrode (ISE)

Turn Around Time: 2-3 days

Reference Range:

Use: Evaluation of renal function, electrolyte balance.

The CPT codes provided are based on AMA guidelines and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.